

AI-Rad Companion Organs RT

Automation in Contouring with AI state-of-the-art algorithms

Ekta Jhala

Sr Manager, SW Product Marketing, EMEA

Varian Medical Systems

General Treatment Planning Fair Balance Safety Statement

Intended Use Summary

Varian Medical Systems' linear accelerators are intended to provide stereotactic radiosurgery and precision radiotherapy for lesions, tumors, and conditions anywhere in the body where radiation treatment is indicated.

Important Safety Information


Radiation treatments may cause side effects that can vary depending on the part of the body being treated. The most frequent ones are typically temporary and may include, but are not limited to, irritation to the respiratory, digestive, urinary or reproductive systems, fatigue, nausea, skin irritation, and hair loss. In some patients, they can be severe. Treatment sessions may vary in complexity and time. Radiation treatment is not appropriate for all cancers.

Medical Advice Disclaimer

Varian as a medical device manufacturer cannot and does not recommend specific treatment approaches. Individual treatment results may vary.

Not all products or features available for sale in all markets

Cancer Affects Everyone



24.6 Million

New Cancer Cases by 2030,
>75% increase from 2010



\$2 Trillion

Global Economic Burden in
2010



50 - 60%

Cancer Patients Requiring RT



56%

Cancers Diagnosed In
High Income Countries



10%

Cancer Patients That Have
Access To RT in Low Income
Countries



> 20,000*

Linacs needed by 2035; greatest
need in low + mid-income
countries


We must improve access to cancer care

By 2035, we will need an additional 150,000 skilled RO clinicians

What is Needed


Radiation Oncology Centers	7,700	3,200	10,900
Linear Accelerators	13,100	21,800*	21,800
Radiation Oncologists	23,200	22,300	45,500
Medical Physicists	10,000	29,300	39,300
Radiation Technologists	33,300	96,900	130,200

GAP



**People, AI and data are key
to closing the gap**

- ✓ Quality
- ✓ Efficiency
- ✓ Simplicity
- ✓ Automation



Source: Expanding global access to radiotherapy. Lancet Oncol. Vol 16, Sept. 2015

* 8,700 new machines + 13,100 replacements = 21,800 machines needed


Challenges with manual contouring



Personnel costs contribute up to 30% of Radiation Therapy Treatment course/patient ¹



Contouring errors are high risk part of the Radiation Oncology Procedure ²



OAR Delineation is still one of the largest sources of variability ³


¹ IAEA, *Radiotherapy in Cancer Care: Facing the Global Challenge*, 2017

² Radiation Oncology Incident Learning System, *Aggregate Report Patient Safety Work Product*, Q4, 2017

³ van der Veen J.a,Gulyban A.b,Willems S.c,Maes F.c,Nuyts S. *Interobserver variability in organ at risk delineation in head and neck cancer*, 2021

Blinded evaluation of auto contouring at Universitätsklinikum Erlangen

Blinded physician rating (%)



- *The case evaluation was conducted with Organs RT on syngo.via RT Image Suite.
- The feedback and the results are from the collaboration performed at UKER
- The statements by Siemens Healthineers' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.

“Current **state-of-the art** enables **OAR auto segmentations** that are **on par with human experts**”

“Important **prerequisite to automate and accelerate RT planning workflow by multiple orders of magnitude (Adaptive Radiotherapy/plan of the day)**”



Dr. Florian Putz


Physicist, Radiation Oncology
Universitätsklinikum Erlangen,
Germany

Study details:

Clinical evaluation of 50 CT datasets for 5 sites with each 10 cases.

3 RT physicians (one Senior Physician, two Physicians) rated auto contouring solutions & peers (each other).


Blinded evaluation of auto contouring at Universitätsklinikum Erlangen (continued)



- *The case evaluation was conducted with Organs RT on *syngo.via* RT Image Suite.
- The feedback and the results are from the collaboration performed at UKER.
- The statements by Siemens Healthineers' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.

Evaluation of autocontouring at Centre de Cancérologie du Grand Montpellier

Blinded physician rating (%)



- 1) Head & Neck (6 datasets), Thorax (25 datasets), Abdomen (2 datasets), Male Pelvis (17 datasets) and Female Pelvis (5 datasets).
- The feedback and the results are from the collaboration performed at CCGM.
- The statements by Siemens Healthineers' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.

77% of contours are **clinically usable** without requiring any edits.

95% of the contours can be used with or without minor edits.


Study details:

55 CT datasets¹ of contours automatically generated by AI-Rad Companion Organs RT


Contour quality evaluated by 2 board certified clinicians (a physician and a physicist)

AI-powered auto contouring deployment


A solution for every customer



AI-powered
OAR auto-segmentation



Fully automated, streamlined workflow



Efficiency

Designed to Achieve

Consistency

Accessibility


Security and Data Privacy



We develop solutions according to the **“Privacy by design and by default”** principle to ensure built-in data protection



We keep your data protected with **centrally managed service operations** for up-to-date software



We empower you to protect your data by **setting your data preferences accordingly**

Important for IT customers



EuroPriSe¹ seal (GDPR)



ISO27001² certified




HIPAA compliance support

¹ teamplay with its core apps Dose, Usage, Protocols, and Images/Images Research has been awarded by EuroPriSe (<https://www.european-privacy-seal.eu/EPS-en/siemens-healthcare-teamplay>)

² The teamplay performance management applications & platform are developed and operated in accordance with an Information Security Management System that was certified for ISO 27001

Value of AI-based autocontouring solutions




Deep learning-trained OAR contouring designed to increase efficiency and consistency

- Courtesy of University Hospital Erlangen, Germany (left image). Eclipse version 17.0 was used to display the auto contouring results (right image).
- Courtesy of Radiologische Allianz Hamburg, Germany (Rendering). Rendering is based on research results that are not commercially available. Future availability cannot be guaranteed.
- The term autocontouring in the context of this presentation means automated contouring of organs-at-risk structures.





Improved workflow efficiency which frees up resources for other planning tasks



Seamless integration into the daily treatment planning workflow

OAR auto contouring supported anatomical sites



Courtesy of University Hospital Aarhus, Denmark (Head&Neck); Courtesy of University Hospital Erlangen, Germany (Abdomen, 3D Rendering);
Courtesy of Leopoldina Krankenhaus Schweinfurt, Germany (Thorax, Pelvis)

OAR auto contouring supported anatomical sites



Deep-learning based autocontouring

Head & Neck

- Brain
- Brainstem
- Eye globe (L/R)
- Lens (L/R)
- Optic nerve (L/R)
- Optic chiasm
- Parotid gland (L/R)
- Submandibular gland (L/R)
- Oral cavity
- Mandible
- Lips
- Larynx
- Glottis
- Supraglottic larynx (L/R)
- Brachial plexus (L/R)
- LN H&N 510K pending

Thorax

- Female breasts (L/R)
- Lung (L/R)
- Lung lobes (RI, RM, RS, LI, LS)
- Whole heart
- Heart chambers (LL, UL, LR, MR)¹
- Endocardium¹
- Aorta
- Individual ribs (24 ribs)
- Sternum

Abdomen

- Liver
- Spleen
- Kidneys (L/R)
- Abdominopelvic cavity

Pelvic area

- Bladder
- Prostate
- Rectum
- Proximal femur (L/R)
- Seminal vesicles
- LN common iliac (L/R)
- LN internal iliac (L/R)
- LN external iliac (L/R)
- LN obturator (L/R)
- LN presacral

Multi-site

- Body contours
- Spinal cord
- Esophagus
- Skeleton

Courtesy of University Hospital Aarhus, Denmark (Head&Neck); Courtesy of University Hospital Erlangen, Germany (Abdomen, 3D Rendering);

Courtesy of Leopoldina Krankenhaus Schweinfurt, Germany (Thorax, Pelvis)

1) Heart chambers and endocardium are model-based and not available with AI-Rad Companion Organs RT.


Advantages of AI-Rad Companion Organs RT

**Powerful platform
combined with long-term
experience in Artificial
Intelligence**



1,100,000,000

**More than 1.1 billion
clinical images
as well as reports, clinical
and genomic data**




**Worldwide super-computing
infrastructure with
700 AI experiments per day**



**More than 650 patent families
related to machine learning and
more than 250 patent families
related to deep learning**

Together we are stronger!

A World Without Fear Of Cancer



* Data on file Varian Medical Systems, Inc 07/2021

Eclipse



Velocity



AI-Rad Companion
Organs RT



Scripting API



RapidPlan



MCO



GPU





ePeer
Review



HyperArc

Elevating Patient Care





A World Without
Fear of Cancer



varian
A Siemens Healthineers Company

(Thank You)

varian

A Siemens Healthineers Company